

# Verifying indices of user-submitted Illumina library pools using Nanopore sequencing and cross-platform demultiplexing

Remi-Andre Olsen<sup>1</sup>, Franziska Bonath<sup>2</sup>, Fanny Taborsak-Lines<sup>2</sup>,  
Sara Sjunnebo<sup>1</sup>, Alfred Kedhammar<sup>1</sup>, Carl Rubin<sup>2</sup>

<sup>1,2</sup>Science for Life Laboratory, National Genomic Infrastructure  
<sup>1</sup>Department of Biochemistry and Biophysics, Stockholm University  
<sup>2</sup>School of Engineering Sciences in Chemistry, KTH Royal Institute of Technology

## 1) Background


The National Genomics Infrastructure (NGI) at SciLifeLab is an accredited core facility offering a wide variety of library prep and sequencing services to researchers.

The highest throughput of the facility consists of **user-prepared Illumina library pools**.

## 2) Problem

The increasing output of Illumina flowcells can improve the time- and cost efficiency of sequencing, provided that the library pools of different users can be sequenced together.

If users do not report the indices of their library pool correctly, this may cause **missing or unexpected indices** as well as **index collisions** between library pools. Thus, the need arises to **confirm the Illumina index identity within and between user-prepared library pools prior to sequencing**.



## 3) Solution

User-prepared Illumina library pools are carried through a Nanopore library prep, where-in each user-prepared pool can be assigned it's own Nanopore barcode.

The barcoded pools are combined and sequenced on a cheap, low-throughput Nanopore flowcell. Each sample can then be demultiplexed on a user project level using the Nanopore barcode and on sample level using the original Illumina index.

Looking at read counts for a given combination of Nanopore barcode and Illumina index allows us to qualitatively and **quantitatively detect index issues**.